




FreeStep™

FreeStep™ is the most complete solution for the assessment of the foot support, biomechanical analysis of lower limb movement and the study of posture.



www.sensormedica.com
info@sensormedica.com



FreeStep™

TECHNICAL FEATURES

SOFTWARE FUNCTIONS

▪ STATISTICS

- Powerful statistics engine
- Data export in format CSV
- Pre-set and free aggregation queries

▪ AUTOMATIC REPORT

- Automatic textual report
- Comparison with normality indices and analytical description of differences than normal
- International texts in 12 different languages
- Instantaneous and modifiable by the operator

▪ COMPARISONS

- Visual and numerical comparison between the various examinations of the same patient.
- Indispensable in the follow up and in the management of pre and post treatment results

COMPATIBLE INTEGRATED DEVICES

Treadmills, Podoscan 2D and Podoscan 3D

WEARABLE COMPATIBLE INTEGRATED DEVICES

Flexinfit, EMG, Moover, Movit Gait

